(I)

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(GB)

- A. Remove the temporary appliance and clean the cavity. Try-in the appliance carefully and proceed with eventual corrections. You can use the special Try-in EnaCem HF bodies available in the same colours of EnaCem HF cement in order to find the right colour: the EnaCem HF Try-in will not cure. After finding the right colour wash away the try-in material.
- B. Clean the surface of the preparation with alcohol or chlorhexidine and sandblast it.
- C. Etch the tooth surface with EnaEtch 37% phosphoric acid for 1 minute.
- D. Apply on the preparation the mixture of **Enabond** and **Enabond Catalyst**. The adhesive should be wiped on the surface with a disposable microbrush. Dry with air to eliminate water and solvent residuals.
- E. Sandblast the internal part of the composite appliance and clean it with alcohol; apply **Enabond** without curing it but dry it carefully.
- F. Apply a small amount of **EnaCem HF** in the colour you have chosen in the internal side of the appliance to be luted, position it on the tooth and condense it mechanically or manually. Remove composite excess and light-cure for 60 seconds from each side of the tooth.

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